



eScreen Electronic Specimen Results Service

Implementation Guide

Version 1.1

March 2019



Contents

About this Guide	3
Audience	3
Overview	3
Background	3
Delivery Details	3
Document Details.....	5
Electronic Results XML Sample	6
Node Definitions	7
1.1 Specimen Result.....	7
1.1 MRO.....	13
1.1 Analyte.....	14
Lookup Codes.....	15
Revision Log.....	19



About this Guide

The eScreen Electronic Specimen Results service describes the implementation of electronic reporting for vendors and clients receiving toxicology results of clinical or toxicology testing performed at U.S. laboratories affiliated with eScreen's network of clinics and national account clients.

Audience

This guide is intended for eScreen integration partners who are developing electronic reporting systems and receive information from eScreen's operations center.

Overview

Background

This document outlines the XML standard, describes available reporting methods and provides some technical details for its implementation. Detailed discussions with eScreen technical staff may be required to schedule and implement these methods.

Delivery Details

The eScreen electronic specimen results service can deliver electronic results, batch or near real-time, via the following methods with the noted security protocols.

Method	Security Protocol	Notes
HTTP	SSL Required	Push only

HTTP

HTTP delivery is currently provided by using XML push technology. A standard set of responses from the integration partner has been developed to facilitate pseudo-transactional delivery. **These responses are required.**

eScreen provides HTTP delivery with partners using these custom XML responses. The content type of the HTTP post is text/xml.

eScreen technical staff can provide code samples and assistance as required for building receiver components.



Response	Description
<eScreenData TransmissionID='xxx' status='S'/>	Successful receipt of transmission. NOTE: xxx is the actual TransmissionID in the Source XML.
<eScreenData TransmissionID='xxx' status='F'/>	Failed receipt of transmission. NOTE: xxx is the actual TransmissionID in the Source XML.

Note: eScreen asks that clients and vendors confirm data has been received to indicate the transmission was successful.
Data should not be validated at this time.

If no successful response is received after up to 30 seconds, eScreen will try to retransmit two more times before a failed notice is reported and the matter is escalated.

Document Details

Document Type	Definition
XML	

Electronic Results XML Sample

The following XML sample demonstrates the Electronic Results XML document that will be delivered.

```
<eScreenData TransmissionID='12345555'>
  <RecordCount>1</RecordCount>
  <Version>1.0</Version>
  <DateReported>08/23/2001</DateReported>
  <TimeReported>08:09:00.1230000</TimeReported>
  <SpecimenResult Disposition="1">
```



```
<eScreenID>12345678</eScreenID>
<ConfirmationNumber></ConfirmationNumber>
<LabLocationCode></LabLocationCode>
<Dilute></Dilute>
<eScreenClientAccount>100002</eScreenClientAccount>
<eScreenClientSubAccount>1234</eScreenClientSubAccount>
<CostCenter></CostCenter>
<LabName>LabOne</LabName>
<LabAccount>1123343</LabAccount>
<ClientName>Mack Trucking</ClientName>
<InternalClientID>AA111</InternalClientID>
<ElectronicClientID>021543468</ElectronicClientID>
<LocationInfo></LocationInfo>
<CollectionSite></CollectionSite>
<CollectionSitePhoneNumber></CollectionSitePhoneNumber>
<CollectionSiteInNetwork>1</CollectionSiteInNetwork>
<ECupCollection>1</ECupCollection>
<AccessionNum></AccessionNum>
<DonorName>Jones, Tom</DonorName>
<SSN>123456789</SSN>
<OtherID>11222</OtherID>
<OtherIDType>D</OtherIDType>
<ExternalDonorID>555555</ExternalDonorID>
<RandomSelectionID></RandomSelectionID>
<DOB>01/01/1920</DOB>
<HomePhone>913-599-0011</HomePhone>
<HomePhoneExt>111</HomePhoneExt>
<WorkPhone>913-599-0011</WorkPhone>
<WorkPhoneExt>222</WorkPhoneExt>
<ChainOfCustody>09223432</ChainOfCustody>
<CollectionDate>08/14/2001</CollectionDate>
<CollectionTime>10:34:00.1230000</CollectionTime>
<LabReceivedDate>08/14/2001</LabReceivedDate>
<LabReceivedTime>10:34:00.9890000</LabReceivedTime>
<LabReportDate>08/21/2001</LabReportDate>
<LabReportTime>09:00:00.123000</LabReportTime>
<VerificationDate>08/22/2001</VerificationDate>
<VerificationTime>10:00:00.343000</VerificationTime>
<ReasonForTest>PE</ReasonForTest>
<SpecimenCollector>Wayne, Bruce</SpecimenCollector>
<SpecimenType>U</SpecimenType>
<PanelID>CD51</PanelID>
<PrintedComments>These are printed comments.</PrintedComments>
<Regulation>N</Regulation>
<RegulationType>FAA</RegulationType>
<COCCConfirmationDate>08/22/2001</COCCConfirmationDate>
<COCCConfirmationOperator></COCCConfirmationOperator>
<MRO>
```



```
<Name>Dr. MRO Test</Name>
<Address>123 Main St.</Address>
<City>Somewhere</City>
<State>KS</State>
<Zip>66210</Zip>
<Phone>9135555555</Phone>
</MRO>
<Analyte Disposition="N">
  <PanelID>CD51</PanelID>
  <AnalyteID>100</AnalyteID>
  <AnalyteName>Marijuana</AnalyteName>
  <SpecimenType>U</SpecimenType>
  <ScreeningCutoffValue></ScreeningCutoffValue>
  <ScreeningCutoffUnit></ScreenincCutoffUnit>
  <ConfirmationCutoffValue></ConfirmationCutoffValue>
  <ConfirmationCutoffUnit></ConfirmationCutoffUnit>
  <Quantity></Quantity>
</Analyte>
<Analyte Disposition="P">
  <PanelID> CD51</PanelID>
  <AnalyteID>200</AnalyteID>
  <AnalyteName>Cocaine</AnalyteName>
  <SpecimenType>U</SpecimenType>
  <ScreeningCutoffValue></ScreeningCutoffValue>
  <ScreeningCutoffUnit></ScreenincCutoffUnit>
  <ConfirmationCutoffValue></ConfirmationCutoffValue>
  <ConfirmationCutoffUnit></ConfirmationCutoffUnit>
  <Quantity></Quantity>
</Analyte>
<Analyte Disposition="N">
  <PanelID>CD51</PanelID>
  <AnalyteID>310</AnalyteID>
  <AnalyteName>Amphetamine</AnalyteName>
  <SpecimenType>U</SpecimenType>
  <ScreeningCutoffValue></ScreeningCutoffValue>
  <ScreeningCutoffUnit></ScreenincCutoffUnit>
  <ConfirmationCutoffValue></ConfirmationCutoffValue>
  <ConfirmationCutoffUnit></ConfirmationCutoffUnit>
  <Quantity></Quantity>
</Analyte>
</SpecimenResult>
</eScreenData>
```



Node Definitions

NOTE: All time fields appear in the following format:

hh:mm:ss.mmmmmmm. Milliseconds may or may not appear.

For example, you may see 10:34:00.1230000 OR 10:34:00 AM

All times are in Central Standard Time (CST).

The following table provides a definition of selected Nodes:

1.1 Specimen Result

Node	Field Length	Value	Valid Characters	Description
eScreenData TransmissionID	36		1234567890 ABCDEF-	
RecordCount		1 or >1	n/a	Real time (HTTP) delivery always has a RecordCount of 1. eScreen sends a single drug test at a time with HTTP. eScreen batch (FTP) delivery can have a RecordCount of >1
Version		1.0	n/a	Current version
DateReported	10	MM/DD/YYYY	n/a	Date this result file was reported.
TimeReported	16	HH:MM:SS	n/a	Time this result file was reported
SpecimenResult Disposition	10	See Disposition Table	n/a	See Disposition Table
eScreenID	10		n/a	The internal eScreen identifier assigned to this result. Same as



Node	Field Length	Value	Valid Characters	Description
				Specimen ID. Will remain constant.
ConfirmationNumber				
LabLocationCode				
Dilute		See Dilute table	n/a	Depending on the value, indicates if the specimen is Dilute or not
eScreenClientAccount	10		n/a	The internal eScreen client account
eScreenClientSubAccount	10		n/a	The internal eScreen client sub account
CostCenter	50		Alpha numeric / special char.	Clients cost center number
LabName	80		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ	Name of the Lab performing the test
LabAccount **	20		n/a	Lab account number
ClientName	40		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234	Name of Client

Node	Field Length	Value	Valid Characters	Description
			56789 &,- '#\()	
InternalClientID	15	VARCHAR	n/a	Your internal Client ID as provided to eScreen during client setup
ElectronicClientID	200	VARCHAR	abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ 0123456789 #@\$%^*(),.- ,	Your ElectronicClientID as provided to eScreen during Client setup
LocationInfo	80	VARCHAR	n/a	Collection location as supplied by the lab
CollectionSite	10	1565 means not available	n/a	Collection location as supplied by eScreen
CollectionSitePhoneNumber	12	XXX-XXX-XXXX		Phone number of Collection Site. If there is no valid phone number the node will be returned empty
CollectionSiteInNetwork		See CollectionSiteInNetwork table	n/a	Depending on value, indicates if collection site is In or Out of Network
ECupCollection		See ECupCollection table	n/a	Depending on value, indicates if collection is an eCup Collection



Node	Field Length	Value	Valid Characters	Description
AccessionNum	15		n/a	Accession number from lab (Currently not in use)
DonorName	40	Last name, First name	abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234 56789	Name of Donor
SSN	15	123456789	n/a	Social Security number of donor
OtherID	20		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234 56789 –	
OtherIDType	1	See OtherIDType table	n/a	See OtherIDType table
ExternalDonorID	20		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234 56789	If ID is passed to eScreen when scheduling an event, ID is passed back with results Note: This value is not always returned to integration partners.



Node	Field Length	Value	Valid Characters	Description
RandomSelectionID	6	123456	n/a	eScreen Batch ID; unique to the time period in which the batch is run. It is not unique to the donor.
DOB	10	MM/DD/YYYY	n/a	Date of Birth of Donor
HomePhone	12	555-444-6666	n/a	Donor's home phone
HomePhoneExt	5	12345	n/a	Donor's home phone Ext.
WorkPhone	12	555-444-6666	n/a	Donor's work phone
WorkPhoneExt	5	12345	n/a	Donor's work phone Ext.
ChainOfCustody	20		abcdefghijklmnopqrstuvwxyz FGHIJKLMNOP OPQRSTUVWXYZ 0123456789'	Custody Control Form number for the test
CollectionDate	10	MM/DD/YYYY	n/a	Date Specimen was collected
CollectionTime	16	HH:MM:SS	n/a	Time Specimen was collected
LabReceivedDate	10	MM/DD/YYYY	n/a	Date lab received specimen
LabReceivedTime	16	HH:MM:SS	n/a	Time lab received specimen
LabReportDate	10	MM/DD/YYYY	n/a	First date lab reported results to eScreen
LabReportTime	16	HH:MM:SS	n/a	First time lab reported results to eScreen



Node	Field Length	Value	Valid Characters	Description
VerificationDate	10	MM/DD/YYYY	n/a	Date Specimen was verified by MRO
VerificationTime	16	HH:MM:SS	n/a	Time Specimen was verified by MRO
ReasonForTest	2	See ReasonForTest table	n/a	Provides reason for test
SpecimenCollector	78			Name of collector if completed via eScreen123
SpecimenType	1	See SpecimenType table	n/a	Type of Specimen collected
PanelID	30		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234 56789	
PrintedComments	500			Any comments MRO would like client to see
Regulation	1	See Regulation table	n/a	See Regulation table
RegulationType		See RegulationType table	n/a	See RegulationType table



Node	Field Length	Value	Valid Characters	Description
COCConfirmationDate	10	MM/DD/YYYY	n/a	When MRO received copy of CCF (Only for DOT)
COCConfirmationOperator	60	Last name, first name	VarChar	Last name, first name

1.2 MRO

Name	80		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234 56789 ,.-'√	Name of MRO
Address	40		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234 56789 ,.#-'/'	Address of MRO
City	32		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234 56789 ,.-'	City of MRO
State	2		n/a	State of MRO



Zip	7		n/a	Zip code of MRO
Phone	10	5554446666	n/a	Phone number of MRO

1.3 Analyte

Analyte Disposition	1	See Disposition (Attribute of Analyte) table	n/a	See Disposition (Attribute of Analyte) table
PanelID	30		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ FGHIJKLMNOP OPQRSTUVWXYZ WXYZ01234 56789	
AnalyteID	10	See AnalyteID (Node Of Analyte) table	n/a	See AnalyteID (Node Of Analyte) table
AnalyteName	20	See AnalyteID (Node Of Analyte) table	n/a	See AnalyteID (Node Of Analyte) table
SpecimenType	1	See SpecimenType table	n/a	Type of Specimen collected
ScreeningCutoffValue				
ScreeningCutoffUnit				
ConfirmationCutoffValue				
ConfirmationCutoffUnit				



Quantity				
----------	--	--	--	--

****LabAccount should not be used to identify clients.**

Lookup Codes

The following tables define eScreen values, which are used instead of text.

Disposition (Attribute of SpecimenResult)

Value	Description
1	Negative
2	Positive
3	Canceled
4	Adulterated
5	Positive UTCD (Unable To Contact Donor)
7	Alcohol Only Neg
8	Alcohol Only Pos
9	Refusal to Test
10	Substituted
11	Policy Violation
12	No Review/No Report
13	Canceled no CCF
16	No Show
17	Quant Only



Dilute (Node of SpecimenResult)

Value	Description
(blank)	Not Provided by Lab / See Printed Comments
0	Not Dilute
1	Dilute

****Dilute node will only be present if the final disposition is a negative or positive.**

OtherIDType (Node of SpecimenResult)

Value	Description
D	Driver's License
P	Passport
M	Military ID
E	Employer ID
O	Other
<i>Empty (No value)</i>	None

ReasonForTest (Node of SpecimenResult)

Value	Description
FU	Follow Up
NS	Not Specified
OT	Other
PA	Post Accident
PE	Pre-employment
RA	Random



RS	Reasonable Suspicion
RD	Return to Duty
PM	Periodic Medical
PR	Promotion
TR	Transfer
AL	All

SpecimenType (Node of SpecimenResult)

Value	Description
U	Urine
B	Breath
H	Hair
O	Other
P	Oral Fluid
L	Blood

Regulation (Node of SpecimenResult)

Value	Description
D	DOT
N	Non DOT



DOT RegulationType (Node of SpecimenResult)

Value	Description
FAA	Federal Aviation Administration
FMCSA	Federal Motor Carrier Safety Administration
FRA	Federal Railroad Administration
FTA	Federal Transit Administration
HHS	US Department of Health & Human Services
NRC	Nuclear Regulatory Commission
PHMSA	Pipeline Hazardous Materials Safety Administration
USCG	United States Coast Guard

Disposition (Attribute of Analyte)

Value	Description
P	Positive
N	Negative
X	Cancelled
<i>Empty (No Value)</i>	Withheld

CollectionSiteInNetwork (Node of SpecimenResult)

Value	Description
0	Out of Network
1	In Network Tier 1
2	In Network Tier 2



3	LabCorp
4	Quest

eCupCollection (Node of SpecimenResult)

Value	Description
1	eCup Collection
0	Not an eCup Collection

AnalyteID (Node of Analyte)

Please contact eScreen Integration team for relevant AnalyteID values.

Revision Log

Value	Descriptor
1.1	- Specimen Result Node - Table 1.1 - Added SpecimenType Node - Added PanelID Node
1.2	Removed mention of FTP, SFTP, and webservice delivery options