



eScreen Electronic Specimen Status Service

Implementation Guide

Version A

March 2018



Contents

About this Guide	3
Audience	3
Overview	3
Background	3
Delivery Details	3
Document Details.....	5
Electronic Specimen Status XML Sample	5
Node Definitions	6
Lookup Codes.....	11



About this Guide

The eScreen Electronic Specimen Status service describes the implementation of electronic status reporting for vendors and clients receiving toxicology services of clinical or toxicology testing performed at U.S. laboratories affiliated with eScreen's network of clinics and national account clients.

Audience

This guide is intended for eScreen integration partners who are developing electronic reporting systems and receive information from eScreen's operations center.

Overview

Background

This document outlines the XML standard, describes available reporting methods and provides some technical details for its implementation. Detailed discussions with eScreen technical staff may be required to schedule and implement these methods.

Delivery Details

The eScreen electronic specimen status service can deliver electronic status updates, batch or near real-time, via the following methods with the noted security protocols.

Method	Security Protocol	Notes
HTTP	SSL Required	Push only

HTTP

HTTP delivery is currently provided by using XML push technology. A standard set of responses from the integration partner has been developed to facilitate pseudo-transactional delivery. **These responses are required.**

eScreen provides HTTP delivery with partners using these custom XML responses. The content type of the HTTP post is text/xml.

eScreen technical staff can provide code samples and assistance as required for building receiver components.



Response	Description
<eScreenData TransmissionID='xxx' status='S'/>	Successful receipt of transmission. NOTE: xxx is the actual TransmissionID in the Source XML.
<eScreenData TransmissionID='xxx' status='F'/>	Failed receipt of transmission. NOTE: xxx is the actual TransmissionID in the Source XML.

Note: eScreen asks that clients and vendors confirm data has been received to indicate the transmission was successful.
Data should not be validated at this time.

If no successful response is received after up to 30 seconds, eScreen will try to retransmit two more times before a failed notice is reported and the matter is escalated.

Document Details

Document Type	Definition
XML	

Electronic Specimen Status XML Sample

The following XML sample demonstrates the Electronic Status XML document that will be delivered.

```
<eScreenData TransmissionID="55eb6020-cc2f-476f-afee-981243f628e0">  
  <Version>1.0</Version>  
  <DateReported>2/1/2017</DateReported>  
  <TimeReported>01:55:38 PM</TimeReported>  
  <SpecimenStatus Status="5">  
    <eScreenID />
```



Abbott

```
<ConfirmationNumber>AI18544185NH</ConfirmationNumber>
<eScreenClientAccount>126234</eScreenClientAccount>
<eScreenClientSubAccount>6</eScreenClientSubAccount>
<LabName>Evidential Breath Test</LabName>
<LabAccount />
<ClientName>DrugScan - Human Resources</ClientName>
<InternalClientID />
<ElectronicClientID>DSI|DRUGSCAN</ElectronicClientID>
<LocationInfo />
<CollectionSite>36183</CollectionSite>
<AccessionNum />
<SpecimenType>B</SpecimenType>
<DonorName>Jones, Jane</DonorName>
<SSN>656565656</SSN>
<OtherID>65656565</OtherID>
<OtherIDType>D</OtherIDType>
<ExternalDonorID />
<DOB>1/1/1965</DOB>
<HomePhone />
<HomePhoneExt />
<WorkPhone />
<WorkPhoneExt />
<ChainOfCustody />
<CollectionDate />
<CollectionTime />
<LabReceivedDate />
<LabReceivedTime />
<LabReportDate />
<LabReportTime />
<VerificationDate />
<VerificationTime />
<ReasonForTest>PE</ReasonForTest>
<SpecimenCollector>Wayne, Bruce</SpecimenCollector>
<PrintedComments />
<Regulation>N</Regulation>
<RegulationType />
  </SpecimenStatus>
<RecordCount>1</RecordCount>
</eScreenData>
```

Node Definitions

NOTE: All time fields appear in the following format:



hh:mm:ss.mmmmmmm. Milliseconds may or may not appear.

For example, you may see 10:34:00.1230000 OR 10:34:00

All times are in Central Standard Time (CST).

The following table provides a definition of selected Nodes:

Node	Field Length	Value	Valid Characters	Description
eScreenData TransmissionID	36		12345678 90ABCDEF -	
RecordCount		1 or >1	n/a	Real time (HTTP) delivery always has a RecordCount of 1. eScreen sends a single drug test at a time with HTTP. eScreen batch (FTP) delivery can have a RecordCount of >1
Version		1.0	n/a	Current version
DateReported	10	MM/DD/YYYY	n/a	Date this status was reported
TimeReported	16	HH:MM:SS	n/a	Time this status was reported
SpecimenStatus Disposition	10	See Disposition Table	n/a	See Disposition Table
eScreenID	10		n/a	The internal eScreen identifier assigned to this status. Same as Specimen ID. Will remain constant.
ConfirmationNumber				Number to identify specific event
eScreenClientAccount	10		n/a	The internal eScreen client account



Node	Field Length	Value	Valid Characters	Description
eScreenClientSubAccount	10		n/a	The internal eScreen client sub account
CostCenter	50		Alpha numeric / special char.	Clients cost center number
LabName	80		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ Z	Name of the Lab performing the test
LabAccount **	20		n/a	Lab account number
ClientName	40		abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ Z0123456 789 &,- '#\()	Name of Client
InternalClientID	15	VARCHAR	n/a	Your internal Client ID created by eScreen during client setup
ElectronicClientID	200	VARCHAR	abcdefghijklmnopqrstuvwxyz ABCDEFGHIJKLMNOPQRSTUVWXYZ LMNOPQR	Your ElectronicClientID as provided to eScreen during Client setup



Node	Field Length	Value	Valid Characters	Description
			STUVWXY Z 01234567 89 #@\$%^*() ,.-'	
LocationInfo	80	VARCHAR	n/a	Collection location as supplied by the lab
CollectionSite	10	1565 means not available	n/a	Collection Site ID as supplied by eScreen
AccessionNum	15		n/a	Accession number from lab <i>(Currently not in use)</i>
Specimen Type				Type of Specimen tested (see Specimen Type table in Lookup Codes section)
Regulation Type				The DOT regulation type (if known)
DonorName	40	Last name, First name	abcdefghij klmnopqrs tuvwxyzAB CDEFGHIJK LMNOPQR STUVWXY Z0123456 789	Name of Donor
SSN	15	123456789	n/a	Social Security number of donor
OtherID	20		abcdefghij klmnopqrs	

Node	Field Length	Value	Valid Characters	Description
			tuvwxyzAB CDEFGHIJK LMNOPQR STUVWXY Z0123456 789 –'	
OtherIDType	1	See OtherIDType table	n/a	See OtherIDType table
ExternalDonorID	20		abcdefghijkl mnopqrs tuvwxyzAB CDEFGHIJK LMNOPQR STUVWXY Z0123456 789	If ID is passed to eScreen when scheduling an event, ID is passed back with status Note: This value is not always returned to integration partners in cases where the ePassport is not used by the donor.
DOB	10	MM/DD/YYYY	n/a	Date of Birth of Donor
HomePhone	12	555-444-6666	n/a	Donor's home phone
HomePhoneExt	5	12345	n/a	Donor's home phone Ext.
WorkPhone	12	555-444-6666	n/a	Donor's work phone
WorkPhoneExt	5	12345	n/a	Donor's work phone Ext.
ChainOfCustody	20		abcdefghijkl mnopqrs tuvwxyzAB CDEFGHIJK LMNOPQR STUVWXY	Custody Control Form number for the test



Node	Field Length	Value	Valid Characters	Description
			Z0123456 789'	
CollectionDate	10	MM/DD/YYYY	n/a	Date Specimen was collected
CollectionTime	16	HH:MM:SS	n/a	Time Specimen was collected
LabRecievedDate	10	MM/DD/YYYY	n/a	Date lab received specimen
LabRecievedTime	16	HH:MM:SS	n/a	Time lab received specimen
LabReportDate	10	MM/DD/YYYY	n/a	First date lab reported results to eScreen
LabReportTime	16	HH:MM:SS	n/a	First time lab reported results to eScreen
VerificationDate	10	MM/DD/YYYY	n/a	Date Specimen was verified by MRO
VerificationTime	16	HH:MM:SS	n/a	Time Specimen was verified by MRO
ReasonForTest	2	See ReasonForTest table	n/a	Provides reason for test
SpecimenCollector	78			Name of collector if completed via eScreen123
PrintedComments	500			Any comments MRO would like client to view
Regulation	1	See Regulation table	n/a	See Regulation table



Node	Field Length	Value	Valid Characters	Description
RegulationType		See RegulationType table	n/a	See RegulationType table

****LabAccount should not be used to identify clients.**

Lookup Codes

The following tables define eScreen values, which are used instead of text.

Status (Attribute of SpecimenStatus)

Value	Description
1	Sent to Lab
2	Received at Lab (NOTE: Not supported by all labs)
3	In MRO Review
4	Pending CCF (DOT Tests Only)
5	Expired (Scheduled Event has Expired)

SpecimenType (Node of SpecimenStatus)

Value	Description
U	Urine
B	Breath
H	Hair
O	Other
P	Oral Fluid
L	Blood



OtherIDType (Node of SpecimenStatus)

Value	Description
D	Driver's License
P	Passport
M	Military ID
E	Employer ID
O	Other
<i>Empty (No value)</i>	None

ReasonForTest (Node of SpecimenStatus)

Value	Description
FU	Follow Up
NS	Not Specified
OT	Other
PA	Post-Accident
PE	Pre-employment
RA	Random
RS	Reasonable Suspicion
RD	Return to Duty
PM	Periodic Medical
PR	Promotion
TR	Transfer



AL	All
----	-----

Regulation (Node of SpecimenStatus)

Value	Description
D	DOT
N	Non-DOT

RegulationType (Node of SpecimenStatus)

Value	Description
FAA	Federal Aviation Administration
FMCSA	Federal Motor Carrier Safety Administration
FRA	Federal Railroad Administration
FTA	Federal Transit Administration
HHS	US Department of Health & Human Services
NRC	Nuclear Regulatory Commission
PHMSA	Pipeline Hazardous Materials Safety Administration
USCG	United States Coast Guard